DIABASS® SecureSend


How to send data from mySugr App



Step 1


Open mySugr App and choose menu REPORTS, select FILE FORMAT "CSV"

Then click on SEND



Step 2

In the appearing list, please choose the **DIABASS SecureSend** app.



DIABASS® SecureSend


How to send data from mySugr App



Step 3

The DIABASS SecureSend App will now start automatically and display the file to be sent.

Please confirm with click on NEXT



Click **SCAN BARCODE** and scan the barcode from the info sheet you received from your practice. The recipient's information will be displayed.

Please confirm with click on **NEXT**

Step 5

In the next step, you can enter an (optional) message for the recipient.

Please click on **NEXT** to proceed.

DIABASS[®] SecureSend

How to send data from mySugr App



Step 6

After some instants, the data are now highly securely encrypted (AES-256) and ready for transmission.

You can now choose how to send the data.

Sending via Email App:

An email message is generated in your default email app. You need to send this E-Mail manually from this email app.

Sending via Server:

The encrypted data will be posted to the *DIABASS SecureSend* server. From there, the data will be forwarded to the recipient via email, using a neutral sender adress. As soon as completed, the data will be deleted from *DIABASS SecureSend* server.

Please click on **NEXT** to proceed.

Proce Pallay	DIABASS
Choose da	ta transfer
Please select transfer m	ethod
Sending via the ema An email is generate you can then send in	d in the mail app which
Sending via server The encrypted data i SecureSend server.	s sent via the DIABASS
Cancel	Next